



**Akzenta**  
**ONE®**

# TOP MASK ULTRA SAFE

## DISPOSABLE SURGICAL 3-LAYERS FACE MASK WITH ELASTIC BAND

### CONFIGURATION OVERVIEW

- Designed for medical use
- Soft and comfortable
- Water repellent
- Fiberglass-free
- Latex-free
- Good performance of anti micro-dust, anti-spray, anti-virus
- Designed for medical use

- Conformity to standards Type-IIR, Bacterial Filtration efficiency (BFE) >99%, differential pressure <60

### TECHNICAL SPECIFICATIONS

- 3 layers (2 polypropylene layers) + 1 Melt blown layer (middle)
- Raw Material
- Outer Layer 35 g/m<sup>2</sup>, PP, Color
- Middle Layer 25 g/m<sup>2</sup>, Purpose-made Melt blown, White
- Inner Layer 25 g/m<sup>2</sup>, PP, nonwoven material, White
- 2 Holding Elastics urethane elastic fiber
- 1 Nose Barret: polyethylene with iron wire
- Force at break 6 N

### QUALITY ASSURANCE

- Tests performed in accordance with EN 14683:2019

### PRODUCT VARIATIONS

Art-No.	FAMAA24-AKZ-WH-K02	FAMAA24-AKZ-HO-K02	FAMAA24-AKZ-PI-K02	FAMAA24-AKZ-GR-K02	FAMAA24-AKZ-FG-K02	FAMAA24-AKZ-BL-K02	FAMAA24-AKZ-MB-K02	FAMAA24-AKZ-LU-K02
Color	Star White	HotOrange	Pink Panther	EmeraldGreen	Fresh Green	Blue	Med Blue	Laguna Blue

Dimensions:	Thickness:	Weight:
17.5 x 9.5 cm	0.6mm	3.05g

## PACKGING

Pieces/dispenser 50 Dispenser/carton 40

Carton size 52x38x34 cm Carton weight 8.5 Kg

Cartons/layer 4 Layers/pallet 6

Cartons/pallet 24

Pallet size 80x120x219 cm Pallet weight 224 Kg

## FEATURES AND INSTRUCTIONS

### Mask Features:

1. Outer Layer:  
Water-repellant, to prevent fluids from entering the mask
2. Middle Layer:  
Being high level static filtration non-woven, effectively preventing micro dust and bacteria with Bacteria Filtration Efficiency (BFE) at 99.55%, Particle Filtration Efficiency (PFE)>99% at 0.1 +10m, tested by Nelson Lab.
3. Inner Layer:  
Soft and comfortable PP non-woven
4. Total Mask:  
Be good at anti-spray, blood penetration resistance at 120mm Hg.

### Instructions:

1. Take and face the outer-surface of the mask.
2. While facing the outer-surface of the mask away from you, loop the ear loops comfortably around your ears.
3. Adjust the noseband to fit nicely on the bridge of your nose, conforming to the contour of your face.
4. Now holding the noseband in place, pull the bottom end of the mask under your chin to fully cover your mouth, your nose, as well as your chin.